
SABDTM Coaxial Kit us biopsy®





Helping Answer Life,s Most Important Questions®



Features	Benefits
A precise match between the SAB	The sample notch is always fully ex-
and the introducer allow exact	posed whether the Trocar or Bevel
placement of the tip and sample	tip is selected.
notch of the needle into the lessio	n.
The hub of the introducer comes	to All Introducer are equipped with a
complete stop against the nose of	threaded luer cap to ensure the sty-
the SABD.	let does not backup as dense tissue
	is entered.
The sample notch is always fully	
exposed whether the 10mm notch	n
or 15 mm notch is selected.	

Cs Kit Product No.	Description	Per Box
SABD-1409-15	SABD 14gx 09cm,20mmT bevel tip 13gx07cm"S" style Introducer	5
SABD-1409-15-T	SABD 14gx 09cm,20mmT trocer tip 13gx07cm"S" style Introducer	5
SABD-1415-15-T	SABD 14gx15cm,20mmT trocer tip 13gx13cm"S" style Introducer	5
SABD-1609-15	SABD 16gx 09cm,20mmT bevel tip 15gx07cm"S" style Introducer	5
SABD-1809-15	SABD 16gx 09cm,20mmT bevel tip 17gx07cm"S" style Introducer	5
SABD-1809-15-T	SABD 18gx 09cm,20mmT trocer tip 17gx07cm"S" style Introducer	5
SABD-1815-15	SABD 18gx 15cm,20mmT bevel tip 17gx13cm"S" style Introducer	5
SABD-1815-15-T	SABD 18gx 15cm,20mmT trocer tip 17gx13cm"S" style Introducer	5
SABD-1820-15	SABD 18gx 20cm,20mmT bevel tip 17gx18cm"S" style Introducer	5
SABD-1820-15-T	SABD 18gx 20cm,20mmT trocer tip 17gx18cm"S" style Introducer	5
SABD-2009-15-T	SABD 20gx09cm,20mmT trocer tip 18gx07cm"S" style Introducer	5
SABD-2015-15	SABD 20gx 15cm,20mmT bevel tip 18gx13cm"S" style Introducer	5
SABD-2020-15-T	SABD 20gx 20cm,20mmT trocer tip 18gx18cm"S" style Introducer	5

The Products, photographs, images and words belong to their respective companies.